



Certificate of Analysis

Product name: CITRATE SYNTHASE (EC 2.3.3.1)

Product code: # CS

Batch: 024.00321

Mfg Date: March 2021

Tests	Specifications	Results	Comment
Activity (U/ml)	> 400	476.5	Passed
Specific Activity (U/mg of proteins) ¹	> 15	36.9	Passed

Unit Unit definition²

One Unit of Citrate synthase is defined as the amount of enzyme required to produce one μ mole of citric acid from oxaloacetic acid and acetyl-coenzyme A per minute at 37°C, under the following assay conditions:

Tris/HCl buffer, pH 8.0	95.00 mM
KCl	500.00 mM
Oxaloacetic acid	0.60 mM
Acetyl-coenzyme A	0.30 mM
DTNB [5,5'-Dithiobis (2-nitrobenzoic acid)]	0.30 mM

References

¹ Determined as *Ref.*: Bradford M. M., *Analytical Biochemistry*, Vol. 72: pp. 248-254 (1976).

² Kurz et al., *Biochemistry*, Vol. 39: pp. 2283-2296 (2000).

Storage and use conditions

The enzyme is supplied as an ammonium sulphate suspension and should be stored at +2/+8 °C. For assay, this enzyme should be diluted in 100 mM imidazole buffer, pH 8.0 containing 300 mM KCl. **Swirl to mix the enzyme suspension immediately prior to use.**

Exp. Date: March 2024

Quality Assurance
Rosaria Cassese

CPC BIOTECH S.r.l.
Sede Legale e Amministrativa/Registered and account office:
Via L. Galvani, 1
20875 Burago di Molgora (MB) – Italia
Tel. +39 039 6898029
Fax. +39 039 6899774

Partita IVA 00751490962
Codice Fiscale 03447450150
R.E.A.1074257
Registro Imprese di Monza e Brianza
Cap. Soc. € 114.000,00 i.v.
www.cpcbiotech.it