

Certificate of Analysis



Lot status: **APPROVED**

MODEL: 1000ml Screw Cap Bottle
MEDIA: Peptone Water 1 + Tween80 0.5% pH 7.2
BATCH: 23C0802
DATE OF MANUFACTURING: 08 Mar 2023

PRODUCT CODE: 001.1/33PSC28.1L
STORAGE CONDITIONS: at 2-25°C until expiry date
EXPIRY DATE: 08 Mar 2024

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless clear solution. White crystals may appear. Swirl to dissolve.	<i>Conforms</i>
pH specification	7.2±0.2	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>
Media composition	Casein Peptone 1.00g/L + NaCl 8.50g/L + Tween80 0.5% Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	<i>Conforms</i>

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

Start analysis	10 Mar 2023	End analysis	24 Mar 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

GHIARONI CODE: 9538281000