

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>1000ml Screw Cap Bottle</b>	<b>PRODUCT CODE</b>	<b>001.1/33PSC28.1L</b>
<b>MEDIA</b>	Peptone Water 1 + Tween80 0.5% pH 7.2	<b>MEDIA COMPOSITION</b>	Casein Peptone 1.00g/L + NaCl 8.50g/L + Tween80 0.5%
<b>BATCH:</b>	<b>23K2702</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING</b>	<b>27 Nov 2023</b>	<b>EXPIRY DATE</b>	<b>27 Nov 2024</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacture

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless clear solution. White crystals may appear. Swirl to dissolve.	<i>Conforms</i>
pH specification	7.2±0.2	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

<b>Start analysis</b>	<b>28 Nov 2023</b>	<b>End analysis</b>	<b>12 Dec 2023</b>
<b>Quality Control</b>	Elisabetta Peri		
<b>Quality Assurance</b>	Alessia Diego		

**GHIARONI CODE: 9538281000**