

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>900ml Bottle with Pierceable Cap</b>	<b>PRODUCT CODE</b>	<b>001.2/33P.900</b>
<b>MEDIA</b>	Peptone Water 1 + Tween80 0.1% pH 7.0	<b>MEDIA COMPOSITION</b>	Casein peptone 1.00g/L + Sodium Chloride 8.50g/l + Tween80 0.1%
<b>BATCH</b>	<b>24A3002</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING:</b>	<b>30 Jan 2024</b>	<b>EXPIRY DATE</b>	<b>30 Jan 2025</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless clear solution.	Conforms
pH specification	7.0±0.2	Conforms
Filling media volume	900ml	Conforms

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours	No growth	Conforms
22.5±2.5°C aerobic for 336 hours		

Start analysis	31 Jan 2024	End analysis	14 Feb 2024
Quality Control	Elisabetta Peri		
Quality Assurance	Alessia Diego		

**GHIARONI CODE: 953808024**