

# Certificate of Analysis



Lot status: **APPROVED**

**MODEL:** 1000ml Screw Cap Bottle  
**MEDIA:** Phosphate Buffer 1 (0.25M) pH 7.2  
**BATCH:** 23G1102  
**DATE OF MANUFACTURING:** 11 Jul 2023

**PRODUCT CODE:** 002.1/33PSC28.1L  
**STORAGE CONDITIONS:** at 2-25°C until expiry date  
**EXPIRY DATE:** 11 Jul 2024

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless clear solution. White crystals may appear. Swirl to dissolve.	<i>Conforms</i>
pH specification	7.2±0.1	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>
Media composition	Monobasic Potassium Phosphate 34.00g/L Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

Start analysis	12 Jul 2023	End analysis	26 Jul 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	