

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>1000ml Screw Cap Bottle</b>	<b>PRODUCT CODE</b>	<b>002.1/33PSC28.1L</b>
<b>MEDIA</b>	Phosphate Buffer 1 (0.25M) pH 7.2	<b>MEDIA COMPOSITION</b>	Monobasic Potassium Phosphate 34.00g/L
<b>BATCH</b>	<b>24B2102</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING</b>	<b>21 Feb 2024</b>	<b>EXPIRY DATE</b>	<b>21 Feb 2025</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless clear solution. White crystals may appear. Swirl to dissolve.	<i>Conforms</i>
pH specification	7.2±0.1	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>
22.5±2.5°C aerobic for 336 hours		

<b>Start analysis</b>	<b>22 Feb 2024</b>	<b>End analysis</b>	<b>07 Mar 2024</b>
<b>Quality Control</b>	Elisabetta Peri		
<b>Quality Assurance</b>	Alessia Diego		