

# Certificate of Analysis



Lot status: **APPROVED**

|                               |                                   |                           |                                       |
|-------------------------------|-----------------------------------|---------------------------|---------------------------------------|
| <b>MODEL</b>                  | <b>1000ml Screw Cap Bottle</b>    | <b>PRODUCT CODE</b>       | <b>002.2/33PSC28.1L</b>               |
| <b>MEDIA</b>                  | Phosphate Buffer 2 (0.05M) pH 6.0 | <b>MEDIA COMPOSITION</b>  | Monobasic potassium phosphate 6.80g/L |
| <b>BATCH</b>                  | <b>24C0403</b>                    | <b>STORAGE CONDITIONS</b> | at 2-25°C until expiry date           |
| <b>DATE OF MANUFACTURING:</b> | <b>04 Mar 2024</b>                | <b>EXPIRY DATE</b>        | <b>04 Mar 2025</b>                    |



Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

| CHEMICAL/PHYSICAL TESTS | SPECIFICATION  | RESULT          |
|-------------------------|--|-----------------|
| Appearance              | Colourless clear solution. White crystals may appear. Swirl to dissolve. | <i>Conforms</i> |
| pH specification        | 6.0±0.1  | <i>Conforms</i> |
| Filling media volume    | 1000ml   | <i>Conforms</i> |

## MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

| TESTS  | SPECIFICATION | RESULT          |
|--|---------------|-----------------|
| 32.5±2.5°C aerobic for 336 hours<br>22.5±2.5°C aerobic for 336 hours | No growth     | <i>Conforms</i> |

| Start analysis    | 05 Mar 2024     | End analysis  | 19 Mar 2024 |
|-------------------|-----------------|---|-------------|
| Quality Control   | Elisabetta Peri |   |             |
| Quality Assurance | Alessia Diego   |  |             |

**GHIARONI CODE: 953821543**