

Certificate of Analysis



Lot status: **APPROVED**

MODEL	1000ml Screw Cap Bottle	PRODUCT CODE	002.2/33PSC28.1L
MEDIA	Phosphate Buffer 2 (0.05M) pH 6.0	MEDIA COMPOSITION	Monobasic potassium phosphate 6.80g/L
BATCH	24D0802	STORAGE CONDITIONS	at 2-25°C until expiry date
DATE OF MANUFACTURING:	08 Apr 2024	EXPIRY DATE	08 Apr 2025

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless clear solution. White crystals may appear. Swirl to dissolve.	<i>Conforms</i>
pH specification	6.0±0.1	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

Start analysis	09 Apr 2024	End analysis	23 Apr 2024
Quality Control	Eisabetta Peri		
Quality Assurance	Alessia Diego		

GHIARONI CODE: 953821543