

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL:</b> 1000ml Screw Cap Bottle	<b>PRODUCT CODE:</b> 002.3/33PSC28.1L
<b>MEDIA:</b> Phosphate Buffer 3 (0.05M) pH 6.8	<b>STORAGE CONDITIONS:</b> at 2-25°C until expiry date
<b>BATCH:</b> 22K2103	<b>EXPIRY DATE:</b> 21 Nov 2023
<b>DATE OF MANUFACTURING:</b> 21 Nov 2022	

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless, clear to slightly hazy, with or without precipitate.	<i>Conforms</i>
pH specification	6.8±0.1	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>
Media composition	Monobasic Potassium Phosphate 6.80 g/L Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

Start analysis	23 Nov 2022	End analysis	07 Dec 2022
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

**GHIARONI CODE: 953821546**