

Certificate of Analysis



Lot status: **APPROVED**

MODEL:	1000ml Screw Cap Bottle	PRODUCT CODE:	002.5/33PSC28.1L
MEDIA:	Phosphate Buffer 5 (0.05M) pH 8.0	STORAGE CONDITIONS:	at 2-25°C until expiry date
BATCH:	23A0401	EXPIRY DATE:	04 Jan 2024
DATE OF MANUFACTURING:	04 Jan 2023		

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless, clear to slightly hazy, with or without precipitate.	<i>Conforms</i>
pH specification	8.0±0.1	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>
Media composition	Monobasic Potassium Phosphate 6.80 g/L Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	<i>Conforms</i>

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

Start analysis	05 Jan 2023	End analysis	19 Jan 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

GHIARONI CODE: 953821548