

Certificate of Analysis



Lot status: **APPROVED**

| | | | |
|-------------------------------|--------------------------------|---------------------------|--|
| MODEL | 1000ml Screw Cap Bottle | PRODUCT CODE | 002.6/33PSC28.1L |
| MEDIA | Buffer B.1 (0.07M) pH 6.0 | MEDIA COMPOSITION | Dibasic Potassium Phosphate 2.00 g/L + Monobasic Potassium Phosphate 8.00 g/L |
| BATCH | 23L0401 | STORAGE CONDITIONS | at 2-25°C until expiry date |
| DATE OF MANUFACTURING: | 04 Dec 2023 | EXPIRY DATE | 04 Dec 2024 |

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

| CHEMICAL/PHYSICAL TESTS | SPECIFICATION | RESULT |
|-----------------------------|--|-----------------|
| Appearance | Colourless, clear to slightly hazy, with or without precipitate. | <i>Conforms</i> |
| pH specification | 6.0±0.05 | <i>Conforms</i> |
| Filling media volume | 1000ml | <i>Conforms</i> |

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

| TESTS | SPECIFICATION | RESULT |
|--|---------------|-----------------|
| 32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours | No growth | <i>Conforms</i> |

| | | | |
|--------------------------|--------------------|------------------------|--------------------|
| Start analysis | 04 Dec 2023 | End analysis | 18 Dec 2023 |
| Quality Control | Elisabetta Peri | <i>Elisabetta Peri</i> | |
| Quality Assurance | Alessia Diego | <i>Alessia Diego</i> | |

GHIARONI CODE: 953821550