

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>1000ml Screw Cap Bottle</b>	<b>PRODUCT CODE</b>	<b>002.7/33PSC28.1L</b>
<b>MEDIA</b>	Buffer B.3 (0.1M) pH 8.0	<b>MEDIA COMPOSITION</b>	Dibasic Potassium Phosphate 16.73 g/L + Monobasic Potassium Phosphate 0.523 g/L
<b>BATCH</b>	<b>24A2202</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING</b>	<b>22 Jan 2024</b>	<b>EXPIRY DATE</b>	<b>22 Jan 2025</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless, clear to slightly hazy, with or without precipitate.	Conforms
pH specification	8.0±0.1	Conforms
Filling media volume	1000ml	Conforms

## MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	Conforms

Start analysis	23 Jan 2024	End analysis	06 Feb 2024
Quality Control	Elisabetta Peri		
Quality Assurance	Alessia Diego		

**GHIARONI CODE: 953821551**