

Certificate of Analysis



Lot status: **APPROVED**

MODEL	1000ml Screw Cap Bottle	PRODUCT CODE	002.8/33PSC28.1L
MEDIA	Buffer B.4 (0.1M) pH 4.5	MEDIA COMPOSITION	Monobasic Potassium Phosphate 13.61 g/L
BATCH	23K2903	STORAGE CONDITIONS	at 2-25°C until expiry date
DATE OF MANUFACTURING	29 Nov 2023	EXPIRY DATE	29 Nov 2024

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless, clear to slightly hazy, with or without precipitate.	Conforms
pH specification	4.50±0.05	Conforms
Filling media volume	1000ml	Conforms

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	Conforms

Start analisys	01 Dec 2023	End analisys	15 Dec 2023
Quality Control	Elisabetta Peri		
Quality Assurance	Alessia Diego		

GHIARONI CODE: 953821552