

Certificate of Analysis



Lot status: **APPROVED**

| | | | |
|------------------------------|--------------------------------|---------------------------|--|
| MODEL | 1000ml Screw Cap Bottle | PRODUCT CODE | 002.9/33PSC28.1L |
| MEDIA | Buffer B.6 (0.7M) pH 6.0 | MEDIA COMPOSITION | Dibasic Potassium Phosphate 20.00 g/L + Monobasic Potassium Phosphate 80.00 g/L |
| BATCH | 23K3003 | STORAGE CONDITIONS | at 2-25°C until expiry date |
| DATE OF MANUFACTURING | 30 Nov 2023 | EXPIRY DATE | 30 Nov 2024 |

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

| CHEMICAL/PHYSICAL TESTS | SPECIFICATION | RESULT |
|-------------------------|--|-----------------|
| Appearance | Colourless, clear to slightly hazy, with or without precipitate. | <i>Conforms</i> |
| pH specification | 6.00±0.05 | <i>Conforms</i> |
| Filling media volume | 1000ml | <i>Conforms</i> |

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

| TESTS | SPECIFICATION | RESULT |
|--|---------------|-----------------|
| 32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours | No growth | <i>Conforms</i> |

| Start analysis | 01 Dec 2023 | End analysis | 15 Dec 2023 |
|-------------------|-----------------|--------------|------------------------|
| Quality Control | Elisabetta Peri | | <i>Elisabetta Peri</i> |
| Quality Assurance | Alessia Diego | | <i>Alessia Diego</i> |

GHIARONI CODE: 953821553