

Certificate of Analysis



Lot status: **APPROVED**

MODEL	1000ml Screw Cap Bottle	PRODUCT CODE	002.9/33PSC28.1L
MEDIA	Buffer B.6 (0.7M) pH 6.0	MEDIA COMPOSITION	Dibasic Potassium Phosphate 20.00 g/L + Monobasic Potassium Phosphate 80.00 g/L
BATCH	24C1102	STORAGE CONDITIONS	at 2-25°C until expiry date
DATE OF MANUFACTURING	11 Mar 2024	EXPIRY DATE	11 Mar 2025

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless, clear to slightly hazy, with or without precipitate.	<i>Conforms</i>
pH specification	6.00±0.05	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

Start analysis	12 Mar 2024	End analysis	26 Mar 2024
Quality Control	Elisabetta Peri		<i>Elisabetta Peri</i>
Quality Assurance	Alessia Diego		<i>Alessia Diego</i>

GHIARONI CODE: 953821553