

Certificate of Analysis



Lot status: **APPROVED**

MODEL:	200ml Screw Cap Bottle	PRODUCT CODE:	004.1/31PSC28.200
MEDIA:	Medium A-Medium 1 pH 6.6	STORAGE CONDITIONS:	at 2-25°C until expiry date
BATCH:	22L2002	EXPIRY DATE:	20 Dec 2023
DATE OF MANUFACTURING:	20 Dec 2022		

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	<i>Conforms</i>
pH specification	6.6±0.1	<i>Conforms</i>
Filling media volume	200ml	<i>Conforms</i>
Media composition	Peptone 6.00g/L + Yeast extract 3.00g/L + Dextrose 1.00g/L + Casein peptone 4.00g/L + Beef extract 1.50g/L + Agar 15.00g/L	<i>Conforms</i>
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	<i>Conforms</i>

Start analysis	20 Dec 2022	End analysis	02 Jan 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

GHIARONI CODE: 953821515