

# Certificate of Analysis



Lot status: **APPROVED**

|                              |                               |                           |  |
|------------------------------|-------------------------------|---------------------------|--|
| <b>MODEL</b>                 | <b>200ml Screw Cap Bottle</b> | <b>PRODUCT CODE</b>       | <b>004.3/31PSC28.200</b>   |
| <b>MEDIA</b>                 | Medium A (EP) pH 7.9          | <b>MEDIA COMPOSITION</b>  | Peptone 6.00g/L + Yeast extract 3.00g/L + Dextrose 1.00g/L + Casein peptone 4.00g/L + Beef extract 1.50g/L + Agar 15.00g/L |
| <b>BATCH</b>                 | <b>23K0303</b>                | <b>STORAGE CONDITIONS</b> | at 2-25°C until expiry date  |
| <b>DATE OF MANUFACTURING</b> | <b>03 Nov 2023</b>            | <b>EXPIRY DATE</b>        | <b>03 Nov 2024</b>   |

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

| CHEMICAL/PHYSICAL TESTS | SPECIFICATION  | RESULT   |
|-------------------------|--|----------|
| Appearance              | Amber coloured, clear to slightly opalescent gel forms in plates | Conforms |
| pH specification        | 7.9±0.1  | Conforms |
| Filling media volume    | 200ml  | Conforms |

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

| TESTS                            | SPECIFICATION | RESULT   |
|----------------------------------|---------------|----------|
| 32.5±2.5°C aerobic for 120 hours | No growth     | Conforms |
| 22.5±2.5°C aerobic for 120 hours |               |          |

|                          |                    |                     |                    |
|--------------------------|--------------------|---------------------|--------------------|
| <b>Start analysis</b>    | <b>06 Nov 2023</b> | <b>End analysis</b> | <b>13 Nov 2023</b> |
| <b>Quality Control</b>   | Elisabetta Peri    |                     |                    |
| <b>Quality Assurance</b> | Alessia Diego      |                     |                    |

**GIARONI CODE: 953821514**