

Certificate of Analysis



Lot status: **APPROVED**

MODEL:	200ml Screw Cap Bottle	PRODUCT CODE:	004.6/31PSC28.200
MEDIA:	Medium F (EP) – Medium 19 (USP) pH 6.1	STORAGE CONDITIONS:	at 2-25°C until expiry date
BATCH:	22K2503	EXPIRY DATE:	25 Nov 2023
DATE OF MANUFACTURING:	25 Nov 2022		

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	6.1±0.1	Conforms
Filling media volume	200ml	Conforms
Media composition	Peptone 9.40g/L + Yeast extract 4.70g/L + Beef extract 2.40g/L + Dextrose 10.00g/L + Sodium Chloride 10.00g/L + Agar 23.50g/L	Conforms
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	Conforms

Start analysis	25 Nov 2022	End analysis	30 Nov 2022
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

GHIARONI CODE: 953821519