

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>200ml Screw Cap Bottle</b>	<b>PRODUCT CODE</b>	<b>004.8/31PSC28.200</b>
<b>MEDIA</b>	Medium 8 (USP) pH 5.9	<b>MEDIA COMPOSITION</b>	Peptone 6.00g/L + Yeast extract 3.00g/L + Beef extract 1.50g/L + Agar 15.00g/L
<b>BATCH</b>	<b>23K0201</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING</b>	<b>02 Nov 2023</b>	<b>EXPIRY DATE</b>	<b>02 Nov 2024</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	5.9±0.1	Conforms
Filling media volume	200ml	Conforms

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours	No growth	Conforms
22.5±2.5°C aerobic for 120 hours		

Start analisys	03 Nov 2023	End analisys	08 Nov 2023
Quality Control	Elisabetta Peri		
Quality Assurance	Alessia Diego		

**GHIARONI CODE: 953821525**