

Certificate of Analysis



Lot status: **APPROVED**

MODEL	200ml Screw Cap Bottle	PRODUCT CODE	004.9/31PSC28.200
MEDIA	Medium H (EP) pH 7.9	MEDIA COMPOSITION	Peptone 5.00g/L + Beef extract 3.00g/L + Agar 15.00g/L
BATCH	24C0502	STORAGE CONDITIONS	at 2-25°C until expiry date
DATE OF MANUFACTURING:	05 Mar 2024	EXPIRY DATE	05 Mar 2025

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	7.9±0.1	Conforms
Filling media volume	200ml	Conforms

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	Conforms

Start analysis	07 Mar 2024	End analysis	12 Mar 2024
Quality Control	Elisabetta Peri		
Quality Assurance	Alessia Diego		

GHIARONI CODE: 953821521