

# Certificate of Analysis



Lot status: **APPROVED**

**MODEL:** 200ml Screw Cap Bottle  
**MEDIA:** Medium 10 (USP) pH 7.2

**PRODUCT CODE:** 004.A/31PSC28.200

**BATCH:** 23A3101

**STORAGE CONDITIONS:** at 2-25°C until expiry date

**DATE OF MANUFACTURING:** 31 Jan 2023

**EXPIRY DATE:** 31 Jan 2024

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	<i>Conforms</i>
pH specification	7.2±0.1	<i>Conforms</i>
Filling media volume	200ml	<i>Conforms</i>
Media composition	Casein peptone 17.00g/L + Soy peptone 3.00g/L + Sodium chloride 5.00g/L + Dextrose 2.50g/L + Dibasic potassium phosphate 2.50g/L + Agar 12.00g/L + Tween80 1%	<i>Conforms</i>
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	<i>Conforms</i>

Start analysis	01 Feb 2023	End analysis	06 Feb 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

**GHIARONI CODE: 953821526**