

Certificate of Analysis



Lot status: **APPROVED**

MODEL: 200ml Screw Cap Bottle
MEDIA: Medium 10 (USP) pH 7.2

PRODUCT CODE: 004.A/31PSC28.200

BATCH: 23F1301

STORAGE CONDITIONS: at 2-25°C until expiry date

DATE OF MANUFACTURING: 13 Jun 2023

EXPIRY DATE: 13 Jun 2024

| CHEMICAL/PHYSICAL TESTS | SPECIFICATION | RESULT |
|-------------------------|---|-----------------|
| Appearance | Amber coloured, clear to slightly opalescent gel forms in plates | <i>Conforms</i> |
| pH specification | 7.2±0.1 | <i>Conforms</i> |
| Filling media volume | 200ml | <i>Conforms</i> |
| Media composition | Casein peptone 17.00g/L + Soy peptone 3.00g/L + Sodium chloride 5.00g/L + Dextrose 2.50g/L + Dibasic potassium phosphate 2.50g/L + Agar 12.00g/L + Tween80 1% | <i>Conforms</i> |
| | Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer | |

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

| TESTS | SPECIFICATION | RESULT |
|--|---------------|-----------------|
| 32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours | No growth | <i>Conforms</i> |

| | | | |
|-------------------|-------------------|--------------------------|-------------|
| Start analysis | 14 Jun 2023 | End analysis | 19 Jun 2023 |
| Quality Control | Claudia Previtali | <i>Claudia Previtali</i> | |
| Quality Assurance | Rosaria Cassese | <i>Rosaria Cassese</i> | |

GHIARONI CODE: 953821526