

Certificate of Analysis



Lot status: **APPROVED**

MODEL	200ml Screw Cap Bottle	PRODUCT CODE	004.A/31PSC28.200
MEDIA	Medium 10 (USP) pH 7.2	MEDIA COMPOSITION	Casein peptone 17.00g/L + Soy peptone 3.00g/L + Sodium chloride 5.00g/L + Dextrose 2.50g/L + Dibasic potassium phosphate 2.50g/L + Agar 12.00g/L + Tween80 1%
BATCH	24D0401	STORAGE CONDITIONS	at 2-25°C until expiry date
DATE OF MANUFACTURING	04 Apr 2024	EXPIRY DATE	04 Apr 2025

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	7.2±0.1	Conforms
Filling media volume	200ml	Conforms

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours	No growth	Conforms
22.5±2.5°C aerobic for 120 hours		

Start analysis	04 Apr 2024	End analysis	09 Apr 2024
Quality Control	Elisabetta Peri		
Quality Assurance	Alessia Diego		

GHIARONI CODE: 953821526