

Certificate of Analysis



Lot status: **APPROVED**

MODEL: 200ml Screw Cap Bottle
MEDIA: Medium 11 (USP) pH 8.3

PRODUCT CODE: 004.B/31PSC28.200

BATCH: 23A3102

STORAGE CONDITIONS: at 2-25°C until expiry date

DATE OF MANUFACTURING: 31 Jan 2023

EXPIRY DATE: 31 Jan 2024

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	8.3±0.1	Conforms
Filling media volume	200ml	Conforms
Media composition	Peptone 6.00g/L + Yeast extract 3.00g/L + Dextrose 1.00g/L + Casein peptone 4.00g/L + Beef extract 1.50g/L + Agar 15.00g/L	Conforms
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	Conforms

Start analysis	01 Feb 2023	End analysis	06 Feb 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

GHIARONI CODE: 953821527