

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>200ml Screw Cap Bottle</b>	<b>PRODUCT CODE</b>	<b>004.B/31PSC28.200</b>
<b>MEDIA</b>	Medium 11 (USP) pH 8.3	<b>MEDIA COMPOSITION</b>	Peptone 6.00g/L + Yeast extract 3.00g/L + Dextrose 1.00g/L + Casein peptone 4.00g/L + Beef extract 1.50g/L + Agar 15.00g/L
<b>BATCH</b>	<b>24D0301</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING</b>	<b>03 Apr 2024</b>	<b>EXPIRY DATE</b>	<b>03 Apr 2025</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	<i>Conforms</i>
pH specification	8.3±0.1	<i>Conforms</i>
Filling media volume	200ml	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	<i>Conforms</i>

<b>Start analysis</b>	<b>03 Apr 2024</b>	<b>End analysis</b>	<b>08 Apr 2024</b>
<b>Quality Control</b>	Elisabetta Peri		
<b>Quality Assurance</b>	Alessia Diego		

**GHIARONI CODE: 953821527**