

Certificate of Analysis



Lot status: **APPROVED**

MODEL: 200ml Screw Cap Bottle
MEDIA: Medium E (EP) pH 7.90

BATCH: 22K2501

DATE OF MANUFACTURING: 25 Nov 2022

PRODUCT CODE: 004.E/31PSC28.200

STORAGE CONDITIONS: at 2-25°C until expiry date

EXPIRY DATE: 25 Nov 2023

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	7.9±0.2	Conforms
Filling media volume	180 ml <i>Add 20 mL of Medium E Additive buffer 10X in 180 mL of Medium E to a final volume of 200 mL</i>	Conforms
Media composition	Casein peptone 5.00g/L + Meat extract 3.00 g/L + Agar 10.00 g/L Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	Conforms

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	Conforms

Start analysis	25 Nov 2022	End analysis	30 Nov 2022
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

GHIARONI CODE: 953821518