

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>200ml Screw Cap Bottle</b>	<b>PRODUCT CODE</b>	<b>004.E/31PSC28.200</b>
<b>MEDIA</b>	Medium E (EP) pH 7.90	<b>MEDIA COMPOSITION</b>	Casein peptone 5.00g/L + Meat extract 3.00 g/L + Agar 10.00 g/L
<b>BATCH</b>	<b>23K0206</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING</b>	<b>02 Nov 2023</b>	<b>EXPIRY DATE</b>	<b>02 Nov 2024</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
<b>Appearance</b>	Amber coloured, clear to slightly opalescent gel forms in plates	<i>Conforms</i>
<b>pH specification</b>	7.9±0.2	<i>Conforms</i>
<b>Filling media volume</b>	180 ml <i>Add 20 mL of <u>Medium E Additive buffer 10x</u> in 180 mL of <u>Medium E</u> to a final volume of 200 mL</i>	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	<i>Conforms</i>

<b>Start analisys</b>	<b>06 Nov 2023</b>	<b>End analisys</b>	<b>13 Nov 2023</b>
<b>Quality Control</b>	Elisabetta Peri		
<b>Quality Assurance</b>	Alessia Diego		

**GHIARONI CODE: 953821518**