

Certificate of Analysis



Lot status: **APPROVED**

MODEL:	1000ml Bottle with Pierceable Cap	PRODUCT CODE:	006.1/33P.1L
MEDIA:	NaCl Peptone Buffer pH 7.0	STORAGE CONDITIONS:	at 2-25°C until expiry date
BATCH:	23H2301	EXPIRY DATE:	23 Aug 2024
DATE OF MANUFACTURING:	23 Aug 2023		

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless, clear to slightly hazy, with or without precipitate.	<i>Conforms</i>
pH specification	7.0±0.2	<i>Conforms</i>
Filling media volume	1000ml	<i>Conforms</i>
Media composition	Casein peptone 1.00 g/L + Sodium chloride 4.30g/L + Potassium dihydrogen phosphate 3.60g/L + Disodium hydrogen phosphate 7.20g/L Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	<i>Conforms</i>

MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

Start analysis	24 Aug 2023	End analysis	07 Sep 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	