

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>900ml Bottle with Pierceable Cap</b>	<b>PRODUCT CODE</b>	<b>006.1/33P.900</b>
<b>MEDIA</b>	NaCl Peptone Buffer pH 7.0	<b>MEDIA COMPOSITION</b>	Casein peptone 1.00 g/L + Sodium chloride 4.30g/L + Potassium dihydrogen phosphate 3.60g/L + Disodium hydrogen phosphate 7.20g/L
<b>BATCH</b>	<b>23J1903</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING:</b>	<b>19 Oct 2023</b>	<b>EXPIRY DATE</b>	<b>19 Oct 2024</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
<b>Appearance</b>	Colourless, clear to slightly hazy, with or without precipitate.	<i>Conforms</i>
<b>pH specification</b>	7.0±0.2	<i>Conforms</i>
<b>Filling media volume</b>	900ml	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

<b>Start analysis</b>	<b>20 Oct 2023</b>	<b>End analysis</b>	<b>03 Nov 2023</b>
<b>Quality Control</b>	Elisabetta Peri		
<b>Quality Assurance</b>	Alessia Diego		

**GHIARONI CODE: 953808005**