

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>900ml Bottle with Pierceable Cap</b>	<b>PRODUCT CODE</b>	<b>006.1/33P.900</b>
<b>MEDIA</b>	NaCl Peptone Buffer pH 7.0	<b>MEDIA COMPOSITION</b>	Casein peptone 1.00 g/L + Sodium chloride 4.30g/L + Potassium dihydrogen phosphate 3.60g/L + Disodium hydrogen phosphate 7.20g/L
<b>BATCH</b>	<b>24A2302</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING:</b>	<b>23 Jan 2024</b>	<b>EXPIRY DATE</b>	<b>23 Jan 2025</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Colourless, clear to slightly hazy, with or without precipitate.	<i>Conforms</i>
pH specification	7.0±0.2	<i>Conforms</i>
Filling media volume	900ml	<i>Conforms</i>

## MICROBIOLOGICAL TEST CONTROLS

- **Sterility control** according to USP/EP

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	<i>Conforms</i>

<b>Start analysis</b>	<b>24 Jan 2024</b>	<b>End analysis</b>	<b>07 Feb 2024</b>
<b>Quality Control</b>	Elisabetta Peri	<i>Elisabetta Peri</i>	
<b>Quality Assurance</b>	Alessia Diego	<i>Alessia Diego</i>	

**GHIARONI CODE: 953808005**