

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL:</b>	<b>900ml Bottle with Pierceable Cap</b>	<b>PRODUCT CODE:</b>	<b>011/33P.900</b>
<b>MEDIA:</b>	Tryptone Soy Broth	<b>STORAGE CONDITIONS:</b>	at 2-25°C until expiry date
<b>BATCH:</b>	<b>23F2801</b>	<b>EXPIRY DATE:</b>	<b>28 Jun 2025</b>
<b>DATE OF MANUFACTURING:</b>	<b>28 Jun 2023</b>		

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	yellow coloured, clear solution	Conforms
pH specification	7.3±0.2	Conforms
Filling media volume	900ml	Conforms
Media composition	Tryptone Soy Broth	Conforms
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	Conforms

- Fertility – Growth Promotion Test** according to USP/EP

Control strain	Medium inoculation level	Incubation Conditions	Recovery Specifications	Recovery Results	Cultural Response
<i>E. coli</i> ATCC 8739	10-100 viable microorganisms	18-72 h at 32.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>P. aeruginosa</i> ATCC 9027	10-100 viable microorganisms	18-72 h at 32.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>S. aureus</i> ATCC 6538	10-100 viable microorganisms	18-72 h at 32.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>B. subtilis</i> ATCC 6633	10-100 viable microorganisms	18-72 h at 32.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>B. subtilis</i> ATCC 6633	10-100 viable microorganisms	18-72 h at 22.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>C. albicans</i> ATCC 10231	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>A. brasiliensis</i> ATCC 16404	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	Positive growth	Positive growth	Conforms

<b>Start analysis</b>	<b>29 Jun 2023</b>	<b>End analysis</b>	<b>13 Jul 2023</b>
<b>Quality Control</b>	Claudia Previtali	<i>Claudia Previtali</i>	
<b>Quality Assurance</b>	Rosaria Cassese	<i>Rosaria Cassese</i>	

**GHIARONI CODE: 953808004**