

# Certificate of Analysis



Lot status: **APPROVED**

**MODEL:** 1000ml Screw Cap Bottle  
**MEDIA:** Buffered Peptone Medium (EP) pH 7.0  
**BATCH:** 23F2802  
**DATE OF MANUFACTURING:** 28 Jun 2023

**PRODUCT CODE:** 012.1/33PSC28.1L  
**STORAGE CONDITIONS:** at 2-25°C until expiry date  
**EXPIRY DATE:** 28 Jun 2024

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Yellow solution. Crystals may appear. Swirl to dissolve.	Conforms
pH specification	7.0±0.2	Conforms
Filling media volume	1000ml	Conforms
Media composition	Peptone 10.00g/L + Sodium Chloride 5.00g/L + Potassium dihydrogen phosphate 1.50g/L + Disodium hydrogen phosphate 3.50g/L	Conforms
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	Conforms

Start analysis	29 Jun 2023	End analysis	13 Jul 2023
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

**GHIARONI CODE: 953820516**