

Certificate of Analysis



Lot status: **APPROVED**

MODEL:	9ml Test Tubes	PRODUCT CODE:	013.2/40PSC16.09
MEDIA:	Neutralizing Cen w Saponin	STORAGE CONDITIONS:	at 2-25°C until expiry date
BATCH:	22C1101	EXPIRY DATE:	11 Mar 2023
DATE OF MANUFACTURING:	11 Mar 2022		

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Light amber coloured, clear to slightly opalescent	<i>Conforms</i>
pH specification	7.0±0.2	<i>Conforms</i>
Filling media volume	9ml	<i>Conforms</i>
Media composition	Neutralizing Cen w Saponin	<i>Conforms</i>
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours	No growth	<i>Conforms</i>
22.5±2.5°C aerobic for 120 hours		

- Fertility – Growth Promotion Test**

Control strain	Medium inoculation level	Incubation Conditions	Recovery Specifications	Recovery Results	Cultural Response
<i>E. coli</i> ATCC 8739	>100 viable microorganisms	45-60 min at 22.5 ± 2.5°C	Good growth	Good growth	Conforms
<i>S. aureus</i> ATCC 6538	>100 viable microorganisms	45-60 min at 22.5 ± 2.5°C	Good growth	Good growth	Conforms
<i>E. coli</i> ATCC 8739	10-100 viable microorganisms	16-20 h at 22.5 ± 2.5°C	Good growth	Good growth	Conforms

Start analysis	15 Mar 2022	End analysis	21 Mar 2022
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	

GHIARONI CODE: 9538013240