

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL:</b>	<b>800ml Screw Cap Bottle</b>	<b>PRODUCT CODE:</b>	<b>063/33PSC28.800</b>
<b>MEDIA:</b>	Sabouraud Dextrose Agar	<b>MEDIA COMPOSITION</b>	Sabouraud Dextrose Agar
<b>BATCH:</b>	<b>24B2001</b>	<b>STORAGE CONDITIONS:</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING:</b>	<b>20 Feb 2024</b>	<b>EXPIRY DATE:</b>	<b>20 Feb 2025</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	5.6±0.2	Conforms
Filling media volume	800ml	Conforms

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	Conforms

- Fertility – Growth Promotion Test** according to USP/EP

Control strain	Medium inoculation level	Incubation Conditions	Recovery Specifications	Recovery Results	Cultural Response
<i>C. albicans</i> ATCC 10231	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	70%≤R%≤200%	111,6 %	Conforms
<i>A. brasiliensis</i> ATCC 16404	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	70%≤R%≤200%	113,0 %	Conforms

Start analysis	20 Feb 2024	End analysis	26 Feb 2024
Quality Control	Elisabetta Peri		
Quality Assurance	Alessia Diego		

**GHIARONI CODE: 953808001**