

Certificate of Analysis



Lot status: APPROVED

MODEL Petri ø90mm Single Wrapped	PRODUCT CODE 1008/10
MEDIA Yeast Extract Glucose + CAF	MEDIA COMPOSITION Chloramphenicol Glucose Yeast Extract Agar
BATCH 89305	STORAGE CONDITIONS at 2-25°C until expiry date
DATE OF MANUFACTURING: 16 Nov 2023	EXPIRY DATE 16 May 2024

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	6.6±0.2	Conforms
Filling media volume	30ml±2ml	Conforms
Packaging	Single Wrapped (S.W.)	Conforms

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours	No growth	Conforms
22.5±2.5°C aerobic for 120 hours		

- Fertility – Growth Promotion Test** according to USP/EP

Control strain	Medium inoculation level	Incubation Conditions	Recovery Specifications	Recovery Results	Cultural Response
<i>S. cerevisiae</i> ATCC 9763	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	70% ≤ R% ≤ 200%	81,4%	Conforms
<i>C. albicans</i> ATCC 10231	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	70% ≤ R% ≤ 200%	93,0%	Conforms
<i>A. brasiliensis</i> ATCC 16404	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	70% ≤ R% ≤ 200%	97,8%	Conforms
<i>E. coli</i> ATCC 8739	≥100 viable microorganisms	72-120 h at 22.5 ± 2.5°C	Inhibited	Inhibited	Conforms

Start analysis	16 Nov 2023	End analysis	21 Nov 2023
Quality Control	Elisabetta Peri	<i>Elisabetta Peri</i>	
Quality Assurance	Celeste Annovazzi	<i>Celeste Annovazzi</i>	