

Certificate of Analysis



Lot status: **APPROVED**

MODEL:	200ml Screw Cap Bottle	PRODUCT CODE:	1067/31PSC28.200
MEDIA:	SDA+CAF	STORAGE CONDITIONS:	at 2-25°C until expiry date
BATCH:	22L1502	EXPIRY DATE:	15 Dec 2023
DATE OF MANUFACTURING:	15 Dec 2022		

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	Amber coloured, clear to slightly opalescent gel forms in plates	Conforms
pH specification	5.6±0.2	Conforms
Filling media volume	200ml	Conforms
Media composition	Sabouraud Dextrose Agar + 0.005% Chloramphenicol	Conforms
	Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer	

MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 120 hours 22.5±2.5°C aerobic for 120 hours	No growth	Conforms

- Fertility – Growth Promotion Test** according to USP/EP

Control strain	Medium inoculation level	Incubation Conditions	Recovery Specifications	Recovery Results	Cultural Response
<i>C. albicans</i> ATCC 10231	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	70%≤R%≤200%	98,9%	Conforms
<i>A. brasiliensis</i> ATCC 16404	10-100 viable microorganisms	48-120 h at 22.5 ± 2.5°C	70%≤R%≤200%	89,8%	Conforms
<i>E. coli</i> ATCC 8739	≥100 viable microorganisms	72-120 h at 22.5 ± 2.5°C	Inhibited	Inhibited	Conforms

Start analysis	15 Dec 2022	End analysis	23 Dec 2022
Quality Control	Claudia Previtali	<i>Claudia Previtali</i>	
Quality Assurance	Rosaria Cassese	<i>Rosaria Cassese</i>	