

# Certificate of Analysis



Lot status: **APPROVED**

<b>MODEL</b>	<b>900ml Bottle with Pierceable Cap</b>	<b>PRODUCT CODE</b>	<b>ME009/33P.900</b>
<b>MEDIA</b>	Fluid Thioglycollate Medium	<b>MEDIA COMPOSITION</b>	Fluid Thioglycollate Medium
<b>BATCH</b>	<b>23J2301</b>	<b>STORAGE CONDITIONS</b>	at 2-25°C until expiry date
<b>DATE OF MANUFACTURING</b>	<b>23 Oct 2023</b>	<b>EXPIRY DATE</b>	<b>23 Oct 2024</b>

Medium is prepared utilizing raw materials declared to be TSE-BSE FREE by the manufacturer

CHEMICAL/PHYSICAL TESTS	SPECIFICATION	RESULT
Appearance	yellow straw coloured, clear to very slightly opalescent solution with upper 10% or less medium pink on standing.	Conforms
pH specification	7.1±0.2	Conforms
Filling media volume	900ml	Conforms

## MICROBIOLOGICAL TEST CONTROLS

- Sterility control**

TESTS	SPECIFICATION	RESULT
32.5±2.5°C aerobic for 336 hours 22.5±2.5°C aerobic for 336 hours	No growth	Conforms

- Fertility – Growth Promotion Test** according to USP/EP

Control strain	Medium inoculation level	Incubation Conditions	Recovery Specifications	Recovery Results	Cultural Response
<i>P. aeruginosa</i> ATCC 9027	10-100 viable microorganisms	24-72 h at 32.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>C. sporogenes</i> ATCC 11437	10-100 viable microorganisms	24-72 h at 32.5 ± 2.5°C	Positive growth	Positive growth	Conforms
<i>S. aureus</i> ATCC 6538	10-100 viable microorganisms	24-72 h at 32.5 ± 2.5°C	Positive growth	Positive growth	Conforms

<b>Start analysis</b>	<b>25 Oct 2023</b>	<b>End analysis</b>	<b>08 Nov 2023</b>
<b>Quality Control</b>	Elisabetta Peri		
<b>Quality Assurance</b>	Alessia Diego		

**GHIARONI CODE: 9538390010**